RFID In-hospital Solutions for Patient Safety and Resource Savings

Increase patient safety in hospital

\oslash Keep it easy for patient

Patient's RFID wristband can be scanned from even over a blanket, thus enabling identity verification before treatment/care without waking patient up.

RFID Wristband



Prevent patient from leaving without permission

Utilise RFID gates to scan patients who pass through them and prevent patients from leaving hospital facilities without permission.



Save time and costs in daily asset management

Asset Management with RFID


- ✓ Manage incoming and outgoing equipment

Leverage RFID technology to effectively manage a wide variety of inventories, assets and other items more efficiently including: • Medical devices such as implants, etc.

- Medical equipment and other assets such as wheelchairs, pumps, etc.
- Waste disposal



CL4NX PE



For more information, please contact your local SATO office, or visit: satoeurope.com


All information in this leaflet is accurate as of 2020 May. Product specifications are subject to change without notice. Any unauthorized reproduction of the contents of this leaflet, in part or whole, is strictly prohibited. All other software, product or company names are trademarks or registered trademarks of their respective owners



Solutions for Hospitals and Medical Institutions



SATO Enhances Safety and Security in Healthcare Facilities



Hospital Storage/Sterilization Center

✓ Manage internal equipment loans with accuracy

Track and manage internal loans accurately using barcode labels attached to each equipment and device.

Streamline usage and inventory management

Efficient replenishment can be achieved by tracking each device/equipment with removable labels or cards printed with barcode.



CL4NX PI

⊘ Achieve accurate & efficient billing

Streamline billing process and ensure bill accuracy by tracking devices used with their removable labels. Medical adminis

Medical administration information sticker

-	
SATO Hospital F4	089037
Product : XXXXXXX Specification : OP-310 60mm*30m Lot No. : 000310	m

Clinical Laboratory/Pathology



♂ Manage specimens in pathology tests using special labels

Labels that are resistant to special types of solvents are used to enable barcode-based management in pathological tests.

•	Pathology	label
•	Pathology	label



03/10/200X



Labels for aliquoted samples and outsourced tests are printed. These labels can be modified for microplates and petri dishes.



Solution Enhance management of blood products for patient safety

We provide solutions for printing compatibility labels, labels for autologous blood and segment labels. Match information on labels with electronic health records for further authentication/verification.

 Compatibility label 		
Duration Nation Description 12345667 Description Desc		
Name: JOE BLOGGS NH: 5555000 DOB: 03/10/150X Male SATO Hospital Adult Emergency Department Patient: A RhD Positive		
Denation No: 1234567 ID: 5550000		
BLOGGS Sample No: Issue No: 0022446588 JOE Time:		

